



Slide 1



Slide 2



Slide 3



Slide 4



Session disclosure:
Smartbox Assistive Technology has developed, and patented licensed technology trademarked as Grid 3 and Grid for iPad. This presentation will focus exclusively on page-sets within Grid 3 and Grid for iPad and will not include information on other similar or related AAC software.

Smartbox

Slide 5

Agenda

- Core-Based Approach to AAC
- Pragmatic-Based Approach to AAC
- Clinical Implications of Core vs. Pragmatic
- What is Grid?
- Super Core and tour
- Voco Chat and tour

Smartbox

Slide 6

Learning Outcomes

First Learning Objective

- Describe one difference between a core-based approach and pragmatic-based approach

Second Learning Objective

- Compare one key difference between Super Core and Voco Chat

Third Learning Objective

- Identify a client population that would benefit from Voco Chat and Super Core, respectively.

Smartbox

Slide 7

Core Versus Pragmatic Organization in AAC

Smartbox


Slide 8

What is Core-Based Approach to AAC?

Use of a small number of words to make unique utterances

Focus is on building sentences with individual words instead of phrases

Emphasis on teaching grammar and morphology



Smartbox

Slide 9

What is Core?


Common words that make up most of what we say


- 100 words account for 60% of what we say in one day
- 200-400 words account for 80% of what we say in one day

Goal to combine words to make longer and more complex utterances

Typically verbs, adjectives, and pronouns

Let's practice...

help

more

Smartbox

Slide 10

who	what	why	where	when	wait	now	then	again
I / me	to	want	come	see	say	this	that	here
my / mine	am is are	stop	go	put	choose	in	on	off
it	can	like	get	good	feel	a	the	time
you / your	do	need	help	more	finish	and	with	place
she / her	have	open	hurt	different	make	some	not	thing
he / him	will	turn	think	same	play	all	up	abc and phrases

Slide 11

What is Core? (cont.)


Communicate for a multitude of reasons

Communicate in a variety of situations

Core words can be flexible

- Turn (positioning)
- Turn (turn-taking activity)
- Turn (turn page in a book)

turn



Slide 12


Expand on Core with Fringe Words

Fringe words are more specific to a topic, environment, or individual

Typically nouns

Build request with Core + Fringe

- I want car
- Car go
- See that car
- Do you have a car





Slide 13

What is a Pragmatic-Based Approach to AAC?

Focus on the how we use language in interactions, not specific words

Language organized according to why one wants to communicate

Communication development versus language development



Slide 14

The Pragmatics Profile for People Who Use AAC (Martin, Small & Stevens, 2017)


Reasons to communicate <ul style="list-style-type: none">• Gaining attention• Directing attention• Requesting• Rejecting• Naming• Commenting• Giving Information• Asking Questions	Participation in conversation <ul style="list-style-type: none">• Managing conversational breakdown• Requesting clarification• Maintaining a conversation• Giving Partner feedback• Changing the topic• Starting, joining, & ending conversations• Using social greetings and partings• Telling and responding to jokes
--	---




Slide 15

Why choose one over the other?


- Access needs
- Fatigue/efficiency
- How much scaffolding is needed?
- Training/learning barriers
- Goals of AAC user



Slide 16



Slide 17





Slide 18




Slide 19

Super Core and Voco Chat





- Developed for symbol communication
- Designed for AAC users AND their communication partners
- Available in Symbol Stix, Widgit and PCS (US English) as well as other languages
- Can be used with a range of access methods
- For use on Windows and iOS devices




Slide 20

Overview




	Super Core	Voco Chat
Sizes available	Super Core 50 and 30 Super Core Learning 12 and 20	12 cells
Suitable age range	All (Child, Teen & Adult options)	All (Child, Teen & Adult options)
Single word vocabulary	Approx. 2,250	Over 1,500
Pre stored phrases	Approx. 250	Over 375
Language organization	Core vocabulary	Pragmatic
Schematic (activity based)	Yes	Yes
Grammar support	Smart grammar (automatic) Magic Wand (manual)	Supported message building Focus on meaning, not grammar



Slide 21

Introducing






Slide 22

Super Core Key Features

- Core vocabulary selected based on research
- Dynamic columns give access to core alongside activity-specific words
- Support for language and literacy development
- Available in both child and teen/adult versions




Smartbox

[illegible]

Slide 23

Super Core Key Features (cont.)

- Motor planning
 - Buttons with specific or related words will always appear in the same spots
 - Consistency of locations helps with efficiency



The screenshot shows the Super Core software interface. It features a grid of buttons, each with a word and a corresponding icon. The buttons are organized into rows and columns. Two arrows point to the 'see' button in the top row and the 'see' button in the bottom row, illustrating the consistency of word placement.

Slide 24

Super Core Key Features (cont.)

Fitzgerald Key

Child-friendly label	Geometrical form	Colour	Examples
Round vowels	Triangle	Yellow	no, there, there, got, better
Action words	Rectangle	Orange	go, run, climb, cut
Naming words	Square	Green	hockey, not, spider
Describing words	Adaptation	Blue	big, small, special, yellow
Little words	Quadrilateral, Diamond, Configuration (not obvious)	Grey	and, of, because, up
Question words	Interlocking Phenomenon	Purple	who, what, why

Smartbox

Slide 25

[illegible][illegible]

Slide 26



Super Core
Tour

Smartbox

Slide 27

Introducing



Voco Chat



Smartbox

Slide 28

Voco Chat Key Features

- Supported message building
- Vocabulary to support independence
- Easy to use and personalize
- Intuitive

The diagram illustrates the Voco Chat interface layout. At the top is the 'Chat writing area'. Below it, on the left, is the 'My toolbar' containing a star icon. To the right of the toolbar is the 'Message pathways' section, which includes icons for 'Cancel', 'Play', 'Stop', 'Repeat', 'Delete', 'Undo', 'Redo', 'Copy', and 'Paste'. Further right is the 'Functions column' with icons for 'Back', 'Forward', 'Home', 'Search', 'Settings', 'Help', 'Log out', and 'Log in'. Below the toolbar and pathways is the 'Quick communication' section with icons for 'Greeting', 'Request', 'Offer', 'Action', 'Location', 'Object', 'Quantity', and 'Time'. To the right of this is the 'Everyday words' section with icons for 'Person', 'Place', 'Thing', 'Action', 'Location', 'Object', 'Quantity', and 'Time'. The 'Smartbox' logo is in the bottom right corner.

Slide 29

Voco Chat Key Features (cont.)

- Focus on participation (e.g. activity grids)

The screenshot shows an activity grid in Voco Chat. The grid contains icons for various activities: 'Play' (a person playing), 'Eat' (a person eating), 'Sleep' (a person sleeping), 'Walk' (a person walking), 'Run' (a person running), 'Jump' (a person jumping), 'Sit' (a person sitting), 'Stand' (a person standing), 'Lie down' (a person lying down), 'Get up' (a person getting up), 'Go to bed' (a person going to bed), and 'Wake up' (a person waking up). The 'Smartbox' logo is in the bottom right corner.

Slide 30

Voco Chat Key Features (cont.)


- Tools to support communication, mental health, and wellbeing
- Programming built-in. Less steps to personalize

The screenshot shows a tablet displaying the Voco Chat interface. The tablet screen shows a grid of icons for communication and wellbeing, including 'Greeting', 'Request', 'Offer', 'Action', 'Location', 'Object', 'Quantity', and 'Time'. The 'Smartbox' logo is in the bottom right corner.

Slide 31



Slide 32



Slide 33



Slide 34

Dedicated hardware



Grid Pad 10s



Grid Pad 12




Grid Pad 15


Smartbox

Slide 35

Lightweight hardware



Talk Pad



Touch Pad

Smartbox

Slide 36


Stay in touch

Andrea.King@thinksmartbox.com

Melanie.Gyilling@thinksmartbox.com

Smartbox

Slide 37



Slide 38